

Bestell-Nr.	Gerätebezeichnung	Geeignet für Stäube der Staubklasse	Empfohlen beim Anfall von Feinstaub
1.667-289.0	NT 75/2 Tact ² Me Tc		Größte Mengen
1.667-288.0	NT 75/2 Tact ² Me		
1.667-287.0	NT 65/2 Tact ² Tc		
1.667-286.0	NT 65/2 Tact ²		
1.667-237.0	NT 55/2 Tact ² Me		
1.667-294.0	NT 75/1 Tact Me Te H		Große Mengen
1.148-437.0	NT 50/1 Tact Te H		
1.148-237.0	NT 30/1 Tact Te H		
1.148-435.0	NT 50/1 Tact Te M		Große Mengen
1.148-235.0	NT 30/1 Tact Te M		
1.148-411.0	NT 50/1 Tact Te L		Große Mengen
1.148-311.0	NT 40/1 Tact Te L		
1.148-211.0	NT 30/1 Tact Te L		
1.148-201.0	NT 30/1 Tact L		
1.667-292.0	NT 75/2 Ap Me Tc		Kleine bis mittlere Mengen
1.667-291.0	NT 65/2 Ap		
1.184-505.0	NT 35/1 Ap		
1.184-503.0	NT 25/1 Ap		
1.378-610.0	NT 22/1 Ap Te L		
1.378-600.0	NT 22/1 Ap L		
1.667-267.0	NT 80/1 B1 M S	 + explosive Stäube	Mittlere Mengen
1.667-239.0	NT 75/1 Me Ec H Z22	 + explosive Stäube	Mittlere Mengen